
































De benodigde inhalatiekracht verschild per inhalator

Inhalatieweerstand conform Incheck DIAL G16 en P. Haidl et al. / Respiratory Medicine 118 (2016) 65e75

1. Zonda		Handihaler		High 30	
2. Easyhaler				High 30	
3. Turbuhaler		Medium		Med High 50	
4. Nexthaler				Med High 35	
5. Spiromax				Medium 30	
6. Genuair				Medium 40	
7. Novolizer				Medium 50	
8. Elpenhaler				?	
9. Diskus				Med Low 30	
10. Ellipta				Med Low 30	
11. Forspiro				?	
12. Breezehaler				Low 50	
13. Cyclohaler				Low 40	
14. Redihaler				Low 20	
15. Respimat				pMDI 20	
16. Dosisaerosol				pMDI 20	